




Prevent clogging of your tips!

Steinle Reagent Filters for BioTek® Dispense Cassettes (1µl, 5µl and 10µl)

- Available for all 1µl, 5µl and 10µl Cassettes.
- Adaptable stainless steel or Nylon mesh sizes between 30µm-200µm (62µm is standard).
- Fully autoclavable and reusable (121°C, 1 bar, 20 min).
- Easy and quickly installed in place of tubing weight.



Filter for 5µl and 10µl Cassettes

No.	Version	Qty	\$
BT-00FL	Single	1	49,-
BT-00FL5	Pack of 5	5	195,-




Filter for all Adapter cassettes (1µl, 5µl and 10µl)

No.	Version	Qty	\$
BT-00AD	Single	1	49,-
BT-00AD5	Pack of 5	5	195,-

Filter for all 1µl Cassettes

No.	Version	Qty	\$
BT-00MP	Single	1	49,-
BT-00MP5	Pack of 5	5	215,-



Manifolds / Adapter for all cassettes (1µl, 5µl & 10µl)

Steinle Manifold for BioTek® Dispense Cassettes

- Splitting 8 silicone tubings into one larger supply tubing.
- For a prolonged distance between system and reagent.
- Ideal for automated environments.
- Supply tubing available with silicone or **PTFE (Teflon®)** in 40 cm- - 200 cm lengths.
 - Teflon tubing is known for non-stick characteristics (therefore more evenly reagent distribution in the well plate).
- Fully autoclavable and reusable.

No.	Version	Qty	\$
BT-SP8ST	For 5µl & 10µl, Single	1	95,-
BT-SP8ST5	For 5µl & 10µl, Pack of 5	5	395,-
BT-SP8SC	For 1µl, Single	1	95,-
BT-SP8SC5	For 1µl, Pack of 5	5	395,-