

ADVANCED

## Manifolds for MicroFill® Dispenser

fully compatible with the  $\mu$ Fill® dispenser from BioTek® / Agilent®







2025

- Identical Heads
- Angled Tip Heads
- Large Bore Heads


**NEW**

- Angled Tip Heads
- Large Bore Heads
- Manifold for 6- well microplates



## ADVANCED


### Steinle Manifolds for MicroFill™ Dispenser 6, 24, 96 and 384 well microplates



#### Cell-friendly angled dispensing

For sensitive cell-based assays, the angled tips allow gentle dispensing to the side of the well.

This helps keep cell layers intact, even through multiple-step operations.



Fluid stream

Intact cell monolayer

#### Identical and customized versions

All Steinle MicroFill® - Manifolds are fully compatible and autoclavable for the MicroFill® dispensers from BioTek® / Agilent®.

*Please contact us for additional customized versions.*

Cat. No.	Description	Qty.	€/ \$ item
ST-40508S	Manifold, 16-channel for 96- and 384 plates	1	POR
ST-40515S	Manifold, 8-channel for 96-well half area and PCR microplates	1	POR
ST-40518S	Manifold, 8-channel for 24-well microplates	1	POR
ST-60508A	<b>Angled tips</b> - Manifold, 16-channel for 96- and 384-well	1	POR
ST-60515A	<b>Angled tips</b> - Manifold, 8-channel for 96-well half area & PCR	1	POR
ST-80516C	<b>Angled tips</b> - Manifold, 8-channel for 6-well plates	1	POR

All prices are excl. tax and any excise duties. The prices are non-binding while errors are excluded.

The calculation takes place in Euro/USD at the prices valid on the day of the conclusion of the contract.

By placing an order, the buyer recognizes and accepts our general terms and conditions, which can be found at our website at [www.steinle-labtec.com](http://www.steinle-labtec.com)

MicroFill® is a registered trademark from BioTek®  
*Agilent®* .

**Steinle Labtechnology**

Ernst-Reuter-Straße 72  
63110 Rodgau  
Germany

CEO Martin Pohl

Phone +49 (0)6106 62524 04  
Fax +49 (0)6106 62524 05

[mail@steinle-labtec.com](mailto:mail@steinle-labtec.com)  
[www.steinle-labtec.com](http://www.steinle-labtec.com)